Impronta definitiva

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Definitive impression

The procedure for taking the impression with the replica can be summarized in the following steps:

1. Positioning of the impression pick-ups (Fig. 22).
2. The replica was unloaded in order not to interfere with the pick-ups while inserting and taking the impression.
3. The pick-ups were reduced in height in order to record the relationship between the arches.
4. Application of the polyester adhesive on the replica and its pick-ups.
5. The impression was taken with the mouth closed (Figs. 23 and 24).
6. Recording of the transfer arch and relationship between the arches.

When the maxillary impression was taken, the impression for manufacture of the mandibular removable partial denture was also taken.

Then the definitive impressions were used to make the master models which were then mounted in the articulator (Fig. 25).

Manufacture of the master model

The master model for this type of work was manufactured with a removable silicone layer around the area of the implants. The model was cast in Rock plaster (Zhermack, Badia Polesine, Rovigo, Italy) and after grinding it (Fig. 26), a mask of the tissue appearance of the model was made using high-precision silicone Platinum 95 (Zhermack, Badia Polesine, Rovigo, Italy) (Fig. 27).

At this point the model was unloaded in the area around the implants (Fig. 28) and the space created was filled by injection with the appropriate silicone (Gingifast Rigid, Zhermack, Badia Polesine, Rovigo, Italy) (Figs. 29 and 30). This system reproduces the peri-implant soft tissues with great reliability, and also facilitates their positioning and removal during subsequent working since the housing surface is free from undercuts.

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