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Fabrication of the master models and preparation of the stabilized plaques

The models were cast using a fourth-class plaster Elite Stone (Zhermack, Badia Polesine, Rovigo, Italy) mixed in a vacuum with distilled water following the manufacturer's instructions. The stabilized plates were fabricated with Elite SC Tray EVO self-hardening resin (Zhermack, Badia Polesine, Rovigo, Italy) after having eliminated the undercuts from the models with wax (Tenatex Red, Kemdent, Swindon, England) and isolated them (Prothyl Isolator Zhermack, Badia Polesine, Rovigo, Italy). The plates obtained in this way reproduce the support bases of the final prosthesis with maximum accuracy, an essential aspect in facilitating the recording of the centricity. The reshaped wax rims (Pam Wax, Milan, Italy) were then positioned on the base plates (Figs. 23 and 24).

Recording of the centricity

The centricity was recorded as follows:

1. Verification of the length of the anterior sulcus.
2. Check of the maxillary occlusal plane (Figs. 25 to 27). This was orientated according to the Camper plane in the sagittal direction. The discrepancy with the old occlusal plane indicated by the previous prosthesis is quite notable (Fig. 28).
3. Correction of the mandibular sulcus based on phonetics, identification of the DVO.
4. Recording of the RC centricity (Fig. 29) and transfer arch.
5. The choice of anterior teeth was based on the patient's photographs and instructions (Fig. 30).

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