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Planning of the mandibular prosthesis

The mask fabricated on the replica immediately indicated how it was impossible to respect the rules of setting a correct cantilever by lengthening the occlusal table up to the first molar (Fig. 59). It was then decided to fabricate an overdenture prosthesis on a mandibular bar with brackets (OT Cap Micro Rhein'83 Srl, Bologna, Italy) positioned on the occlusal surface. Even though this type of prosthesis increases the overall dimensions of the prosthetic body compared to the fixed solution, it provides excellent stability during function, with the advantages of easy removal for hygiene maintenance and increased extension of the occlusal table.

Mandibular bar manufacture

A wax-up of the structure was made, the volumes being verified using the vestibular and lingual masks and the type of anchorage to be used being selected (Fig. 60). Then the wax-up was scanned for CAM manufacture from which a bar that reflects the initial design in detail with maximum reliability is obtained (Fig. 61).

The housings for the retentive components were threaded on the occlusal surface of the structure. An additional check with the mask of the diagnostic set-up was then used to verify the spaces available for the superstructure and the teeth (Fig. 62).

In this case the superstructure was also made using the technique of casting on a duplicate model just as used for the maxillary prosthesis, and the housings of the Teflon copings were reproduced directly in the plate in order to reduce the overall dimensions (Fig. 63). High-resilience retentive components were chosen (soft Copings OT Cap Micro Rhein'83 Srl, Bologna, Italy) to promote the distribution of the support between the implants and mucosa.

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